


ASSIGNING AND VIEWING SURVEY RESULTS

ASSIGNING THE ASSESSMENT/SURVEY

1. Log into the Admin panel from your web browser (Chrome)
<https://admin.remetrichealth.com>
2. Select the patient, and select edit
3. In the assessment field choose the assessments to show on the tablet
4. Save the patient page.
5. The patient now has the ability to complete the survey on the tablet.



Hannah Novy-Russell - Edit

Login as patient

PIN NUMBER 111111 [reset](#)

FIRST NAME Hannah

LAST NAME Novy-Russell

EMAIL

AFFILIATION [PMD-Lansdale](#)

TABLET ALERTS On

Affiliate Information

Affiliation*

Tablet Alerts

Spiro PD

Tablet

Managers

Primary diagnosis

Secondary diagnosis

Account assessments

TAKING THE SURVEY ON THE TABLET



1. From the tablet dashboard, select My Symptoms

RemetricHealth


2. Select the survey to take if there are more than one. Otherwise tap the correct answer for the question presented.
3. Tap Next to view and answer the next question until complete.
4. The associated score is shown on the tablet dashboard and results available on the patient online profile.



REVIEWING THE SURVEY RESULTS AND RESPONSES


1. Log into the provider dashboard
 - a. Through the admin panel by selecting login as from your provider profile
 - b. Directly logging into the provider portal <https://admin.remetrichealth.com>

2. Select the patient name in the provider dashboard to see page.



RemetricHealth

3. Once in the patient profile, click on the Symptoms tab to view the calendar. Survey scores are shown on the calendar.
4. Click on the score to review the details of the survey.



If you have any questions, please contact customer@remetrichealth.com.